

### Picus parts and materials



On/off button( POM)

ID Window (PC)

USB Hatch (POM)

Start button (PP)

Control wheel (PP/EPDM)

Charge contact/cover (Gold plated brass/POM)

Tip ejector button (POM)

Keyboard (PA)

Display (PC)

Handle (ASA/PC)

Tip ejector collar  
10 $\mu$ l, 120 $\mu$ l, 300 $\mu$ l, 1000 $\mu$ l (PVDF)  
5000 $\mu$ l, 10ml (PP)

Piston  
10 $\mu$ l (SS)  
120 $\mu$ l, 5000 $\mu$ l, 10ml (PPS)  
300 $\mu$ l, 1000 $\mu$ l (PEI)

Safe-Cone standard filter (PE)

Seal  
10 $\mu$ l, 10ml (FKM)  
120 $\mu$ l, 300 $\mu$ l, 1000 $\mu$ l, 5000 $\mu$ l (EPDM)

Tip cone  
10 $\mu$ l, 120 $\mu$ l (PVDF)  
300 $\mu$ l, 1000 $\mu$ l, (PEI)  
5000 $\mu$ l, 10ml (PPS)



Tip cone housing (PA)

Piston  
10 $\mu$ l (SS)  
120 $\mu$ l, 300 $\mu$ l, 1200 $\mu$ l (PEI)

Seal  
10 $\mu$ l (FKM)  
120 $\mu$ l, 300 $\mu$ l, 1200 $\mu$ l (EPDM)

Tip ejector bar  
10 $\mu$ l, 120 $\mu$ l, 300 $\mu$ l (PA)  
1200 $\mu$ l (LCP)

Tip cone  
10 $\mu$ l, 120 $\mu$ l (PVDF)  
300 $\mu$ l (PEI)  
1200 $\mu$ l (PPS)